

PULMONX

PRODUCT OVERVIEW

Zephyr® Endobronchial Valves

Valves are available in three different sizes.

Description	Catalog number
Zephyr 4.0 Endobronchial Valve	EBV-TS-4.0
Zephyr 4.0-LP Short Endobronchial Valve	EBV-TS-4.0-LP
Zephyr 5.5 Endobronchial Valve	EBV-TS-5.5



Zephyr® Delivery Catheters

Description	Catalog number
Zephyr 4.0* Delivery Catheter	EDC-TS-4.0
Zephyr 4.0-J* Delivery Catheter	EDC-TS-4.0-J
Zephyr 5.5 Delivery Catheter	EDC-TS-5.5



* Note: This catheter can be used for all valves with a diameter of 4.0 mm.

Chartis System

Description	Catalog number
Chartis Pulmonary Assessment System - Console	CHR-CO-300
Chartis Catheter	CHR-CA-12.0
Chartis Catheter XL**	CHR-CA-12.0-XL

** Note: the Chartis Catheter XL is identical to the Chartis Catheter with the exception that the strain relief is removed. This allows the catheter to be inserted approximately two inches further into the bronchoscope biopsy valve. The physician should use extra caution while manipulating the catheter in the biopsy valve area due to the lack of a strain relief.



PULMONX

PRODUCT OVERVIEW

Product information for ZEPHYR® Endobronchial Valves, Delivery Catheters and CHARTIS Pulmonary Assessment System



Pulmonx International Sàrl
Rue de la Treille 4
2000 Neuchâtel
Switzerland
Tel: +41 32 475 20 70
Fax: +41 32 475 20 71
info@pulmonx.com

Pulmonx United Kingdom
Tel: 0808 234 03 76
Fax: 0808 234 03 77
uk-info@pulmonx.com

Pulmonx Australia Pty Ltd.
Tel: 1300 785 666
Fax: 1800 785 666
customerservice@pulmonx.com.au



www.pulmonx.com